

Introducing....

Micro Volume Pipettor

- Fully Autoclavable
- UV Resistance
- Body grip made of TPE (Thermo Plastic Elastomer)
 - Minimum transfer of body heat to the pipette resulting in high accuracy even on continuous use
 - Minimum fatigue of hand & fingers
 - Reduced thumb stress
- Easy identification through color coding
- Very attractive looks
- Ergonomic Design ensuring light weight 7 soft plunger movement.
- Click - stop digital system for volume setting.
- Colour codes for easy identification.
- Calibration conform to DIN 12650 & EN-ISO 8655 standards ensuring high accuracy & precision.
- Calibration report enclosed with every pipette.





**Single Channel
Variable Volume Micropipette**

Cat. No.	Volume Range [μl]	Increments [μl]
CV-100	0.5-10	0.1
CV-200	5-50	0.5
CV-300	50-200	1
CV-400	200-1000	5
CV-500	10-100	0.5
CV-600	100-1000	5
CV-700	2-20	0.1
CV-800	20-200	1
CV-900	1000-5000	50
CV-905	500-5000	50
CV-1000	1000-10000	100
CV-1005	500-10000	100

**Single Channel
Fix Volume Micropipette**

Cat. No.	Volume [μl]
CF-5	5
CF-10	10
CF-20	20
CF-25	25
CF-50	50
CF-100	100
CF-200	200
CF-250	250
CF-500	500
CF-1000	1000
CF-2000	2000
CF-5000	5000
CF-10000	10000

8 Channel Micropipette

Cat. No.	Volume Range [μl]	Increments [μl]
CM8-010	0.5 - 10	0.1
CM8-050	5 - 50	0.5
CM8-100	10 - 100	0.5
CM8-200	20 - 200	0.5
CM8-300	30 - 300	1.0

12 Channel Micropipette

Cat. No.	Volume Range [μl]	Increments [μl]
CM12-010	0.5 - 10	0.1
CM12-050	5 - 50	0.5
CM12-100	10 - 100	0.5
CM12-200	20 - 200	0.5
CM12-300	30 - 300	1.0

136 A vieux ain naadja

Algeirs 16000 Algeria

Phone : 023.52.03.91 Fax : 023.52.03.92

Website : www.celabo-dz.com

