

Find us on



Magnus



▼ Theia-fi

MICROSCOPE MASTERCLASS

- Parfocal and centered SP achromatic optics
- Choice of halogen, led illumination & battery backup
- Now available in plan optics also
- Optics with multi-layer coating
- Easy access for lamp replacement from front

Optional Accessories



Phase Contrast
Attachment



Darkfield
Attachment



Magnus Opto Systems India Pvt Ltd
Corporate Office: A-5, Mohan Co-operative Industrial Estate,
Mathura Road, New Delhi - 110044



+91-120-4646900



www.magnusopto.com



SPECIFICATIONS

Magnus Inclined Biological Microscope Model Theia-fi

Item	Specifications	Theia-fi Binocular (SP)	Theia-fi Trinocular (SP)										
Optical System	Semi-Plan (SP) Optics	✓	✓										
Body	Aluminium Die- cast body based on ball bearing and wire guides thereby ensuring smooth and precise manipulation	✓	✓										
Inclined Observation Head	Binocular 45 degree inclined, 360 degrees rotatable	✓											
	Trinocular Head 45 degree inclined, 360 degrees rotatable		✓										
Eyepiece (widefield) For Observation	HWF 10x (FN 18mm) paired eyepiece. The unique optical design of the compensating eyepiece provides relief from eye fatigue and renders color-compensated wide-field images of utmost clarity	✓	✓										
Nosepiece	Quadruple revolving nosepiece based on precision ball-bearing mechanism with positive click stop	✓	✓										
Objectives	<table><tr><td><u>Achromat Objectives</u></td><td><u>N.A.</u></td></tr><tr><td>4x</td><td>0.10</td></tr><tr><td>10x</td><td>0.25</td></tr><tr><td>40x (spring loaded)</td><td>0.65</td></tr><tr><td>100x (oil, spring loaded)</td><td>1.25</td></tr></table> Uniformly centered, interchangeable & parfocal Anti-fungus treated Tropicalized anti fungus treatment ensures image excellence for long periods in conditions favouring to fungus growth	<u>Achromat Objectives</u>	<u>N.A.</u>	4x	0.10	10x	0.25	40x (spring loaded)	0.65	100x (oil, spring loaded)	1.25	✓	✓
<u>Achromat Objectives</u>	<u>N.A.</u>												
4x	0.10												
10x	0.25												
40x (spring loaded)	0.65												
100x (oil, spring loaded)	1.25												
Mechanical Stage	Co-axial low drive mechanical stage (125mm x 145mm) (+/-5mm) with traverse area of 30mm x 76mm (+/-5mm) with single slide holder	✓	✓										
Condenser Holder	Rack and pinion mounted condenser holder	✓	✓										
Condenser	Pre-centered abbe condenser with aperture iris diaphragm (N.A. 1.25) focusable with rack & pinion through 10mm and a continously variable iris diaphragm with a removable blue filter for daylight observation	✓	✓										
Focusing System	Co-axial coarse and fine controls with a focus adjustment and fine adjustment knobs. Coarse focus range 20mm. Fine focus rotation 0.2mm	✓	✓										
Illumination Base With Option	(a) Built in illumination base with pre-centered 6V 20W halogen light source coupled with an efficient collector lens system. Universal power supply 100V-230V AC 50Hz single phase (SMPS)	✓	✓										
	(b) Built in illumination base with pre-centered 6V 20W halogen light source coupled with an efficient collector lens system. Power supply 230V AC 50Hz	✓	✓										
	(c) 3W LED light soure high brightness, longlife (30,000 hrs)	Theia-fi LED (SP)	Theia-fi LED (SP)										
	(d) 1 W LED light soure High brightness, longlife (30,000 hrs). Battery back-up in built NiMH Rechargeable batteries	Theia-fi Freedom (SP)	Theia-fi Freedom (SP)										
● Packed in a corrugated box, with operational manual, dust cover, power card & immersion oil (5ml) ● The reflector mirror, spare bulb & wooden cabinet are optionally available													

